# **New Cooltech Define® Cryoadipolysis Applicators: A Scientific and Comparative Study with Cooltech® Applicators**

---

Eurasian Journal  
of Medical Investigations.  
2019;3(2):85–94  
Gregorio Viera Mármol,  
Jorge Villena

# Key takeaway

- Comsol Multiphysics simulations were used to evaluate the designs of new Cooltech Define applicators, and how changing shape, dimensions and cooling temperature can increase fat loss.
- Variables assessed were: Cooling merit parameter, cooling homogeneity, cooling dynamics of skin and fat, the time required to reach crystallization temperatures (set at 10°C), fat percentage inside the applicators, and time to hypoesthesia.
- Cooltech Define applicators shows better performance compare with Cooltech (the 1st generation):
  - o Quicker and more homogenous cooling.
  - o Faster reach to numb temperature leads to greater patient comfort.
  - o Faster reach to crystallization temperature and longer apoptosis time leads to better result.

**Table 4.** Comparative table of mean times at which crystallization temperatures were reached

Series	Applicator Model	Time for fat to reach crystallization temperature (10.38 °C)
Cooltech®	Straight	48
	Tight	26
	Double	28
Cooltech Define®	Straight	25
	Tight	17
	Double	19


°C: degrees Celsius.

**Table 5.** Comparative table of mean times at which hypoesthesia temperatures were reached

Series	Applicator Model	Time required for 50% of the skin to reach a temperature below 7 °C (min)
Cooltech®	Straight	26
	Tight	9
	Double	8
Cooltech Define®	Straight	3
	Tight	3
	Double	3

°C: degrees Celsius; min: minutes]

Comparison table between cooltech vs cooltech define applicators and the time needed to reach fat crystallization.

Internal temperature distribution in a new cooltech define applicator, showing low-fat temperature and high cooling homogeneity.

# **Safety, Efficacy and Patient Tolerance of the New 360° Cryoadipolysis Device for Multiple-Area Fat Removal**

---

Eurasian Journal of Medical  
Investigations. 2020; 4(1): 35-42

Gregorio Viera Mármol,

Carmen Cano Ochoa,

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Jorge Villena Garcia,

Reyna Vargas Lamas,

Mar Altaba Rosas,

Núria Adell Gómez,

Jackie Law

# Key takeaway

- A study was carried out on 33 patients with localized fat using simultaneous multi-procedure was performed on up to 4 areas with different cryoadipolysis applicators.
- The new generation Cooltech Define has a new 360° homogenous cooling plate design and features 4 applicators for providing up to 4 simultaneous treatments.
- They decrease tissue temperature more quickly than the previous generation cryoadipolysis device and 36% quicker to reach cold hypoesthesia.
- After one cryoadipolysis session, the average skin fold reduced 7.80mm with a maximum reduction was 20mm or 50% reduction.
- The average circumference reduction was 32.3 mm with the highest result achieved in this study is 80 mm.
- Comparing with the previous generation cryoadipolysis device, Cooltech Define achieved greater results for localized fat reduction while simultaneously providing greater patient comfort.

**Table 4.** Mean value of the measurements obtained before and after treatment in study 2: weight, skin fold and perimeter circumference

Variable	Pre-procedure Mean (SD; range) n=69	Post-procedure Mean (SD; range) n=69	Difference Mean (SD; range) n=69	p
Weight (kg)	68.64 (12.75; 47–102)	67.84 (12.81; 45–102)	-0.80 (1.35; (-5)–0.80)	5.99E-06
Skin fold (mm)	41.97 (8.89; 16–55)	34.17 (8.60; 12–55)	-7.80 (3.46; (-20)–(0))	1.62E-28
Perimeter circumference (mm)	935.9 (124.9; 550–1180)	903.6 (123.1; 550–1140)	-32.3 (18.6; (-80)–(0))	2.20E-22

*Iconography of the position using a simultaneous application of 4 applicators in the same patient: 1 Oval Curved applicator on the upper abdomen, 1 Straight applicator on the lower abdomen and an Oval Curved applicator on each side of the flanks.*

*The follow-up at 57 days recorded a 50 mm circumference reduction on the lower abdomen and 45 mm on the flanks; the skin fold reductions were 13 mm (36%) on the lower abdomen, 10 mm (31%) on the left flank and 11 mm (33%) on the right flank.*

*Before and after results using a simultaneous application of 3 applicators: 1 Straight applicator on the central infraumbilical abdomen, 1 Oval applicator on the left outer thigh and 1 Oval applicator on the right outer thigh. Follow-up at 120 days: the circumference reduction on the central infraumbilical abdomen was 75 mm; the skin fold reduction was 20 mm or 50%. On the left and right outer thighs, the circumference measurement reduced by 70 mm and the skin fold by 15 mm or 30%.*

*Before and after results using a simultaneous application of 2 Oval applicators, 1 on the left outer thigh and 1 on the right outer thigh. The follow-up at 77 days presented a 40 mm circumference reduction; the skin fold reductions were 10 mm (left) and 10 mm (right).*

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